A:07/09

VIPI FLASH

VIPI FLASH

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VIPI FLASH

VIPI FLASH

Self-curing denture base acrylic resin (powder and liquid) for repairs, rebasing and additions, according to ISO 1567:1999, Type 2, Class 1. VIPI FLASH resins are especially formulated to obtain easy, fast and precise results.

Instructions:
SHAKE THE POWDER-RECIPIENT BEFORE USE.

Preparation:

Before starting any repair, rebasing or addition make sure that the involved pieces are free from any grease or wax residues, to assure a perfect chemical adhesion between the resins.

We recommend an ambient temperature between 15°C and 30°C for mixture and application.

Dose the necessary amount of VIPI FLASH liquid to the glass or silicone mixing recipient, depending on the quantity being prepared. Add the desired VIPI FLASH resin up to an uniform mixing, without powder or liquid excess; If excess exists, discard the extra material and mix the dough, with a spatula until total homogenization, to eliminate air bubbles. Make sure that the mixture has a viscous fluid consistency.

Do not use the resin after it loses the fluid consistency (5 - 6 minutes after dosage, at 23°C).

Do not touch the repaired surfaces with the fingers.

Polymerization:

To polymerize the resin take it to a pressure polymerizer filled with the necessary amount of water, pressurize at 1.5 bar (20 psi) and complete the polymerization for 20 minutes.

Wait until polymerization is completed to start finishing and polishing.

For repairs, additions and rebasing follow the usual procedures.

This resin may be used for orthodontics but it is better to use ORTO COR resin, especially formulated for that purpose.

The VIPI FLASH shades are highly stable and the resin is highly resistant to wear stress.

Note: VIPI recommends for this product the use of mechanical polishing, only.

Precautions:

• Follow usual procedures for infection-control, waxing-up, articulation and maintenance of the prosthesis, as well as for individual protection;

• For dental professional use only;

• Liquid: Non-toxic very inflammable material (Flash Point: 10°C / 50°F);

• Powder: May form inflammable resin / air mixtures;

• Do not allow unpolymerized material to contact the mouth tissue. The polymerization heat may cause burnings;

• For all acrylic resins, handling in ventilated places, using nitrile-gloves and eye-protection, is highly recommended;

• The product should not be handled during pregnancy, lactancy or by known methacrylate sensitive individuals;

• Store the product in a cool place, protected from direct sunlight;

• Do not allow the products to reach the sewage system, nor discard them with the domestic waste.

Composition:

Powder:	Liquid:
• Polymethylmethacrylate	• Methylmethacrylate
• Benzoyl Peroxide	• DMT
• Biocompatible Pigments	• Inhibitor

Shades:

Powder:
Clear, Pink and Black Pink

Packages:

Powder: Containers of 20 g, 80 g, 225 g, 450 g, 1 kg, 2.250 kg and 5kg.

NOTE: Black Pink shade only in 80g.

Liquid: Bottles of 50 ml, 120 ml, 250 ml, 500 ml, 1 liter and 5 liters (Can).

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Madespa S.A.